

### CERO Setting System out of Operation, no use for > 4 weeks

- **Preparing and Performing System Shutdown:**

1. Open the Top Sliding Cover.
2. Make sure that the CERO is switched off before performing the decontamination procedure.
3. Verify that the plugs of the reader are protected. Cap them if necessary.
4. ***Spraying directly on the unprotected reader plugs may cause damage and may cause an electric short of the device.***
5. Spray the CERO Culture Chamber with 70% EtOH disinfectant spray and wipe using a lint free cloth.
6. Close the top sliding cover.
7. Disconnect the main power cable from the power supply.
8. Switch off the CO2 supply.

**System can be stored for an indefinite amount of time!**

- **System Restart:**

1. Install main power cable
2. Connect the CERO to the low pressure side of a CO2 regulator
3. Set the output pressure of the regulator to  $1.0 \pm 0.2$  bar ( $14.5 \pm 3$  psi).
4. Make sure that the connection is tight and secure.
5. Turn the CERO ON Power Switch on the side of the instrument.
6. Open the top sliding cover.
7. Make sure that the CERO is switched off before performing the decontamination procedure.
8. Verify that the plugs of the reader are protected. Cap them if necessary.
9. ***Spraying directly on the unprotected reader plugs may cause damage and may cause an electric short of the device.***
10. Clean surface by whipping with 70% EtOH
11. Connect the Mains Power Cable to the Mains Power Connector on the rear of the CERO.

**System is ready for use!**

