

HEKin1™

Serum free and animal
component-free medium designed
for HEK293 cells



HIMEDIA®

For Life is Precious

 **HIMEDIA**®
Cell Culture
Enabling Breakthroughs

Product Portfolio

Culturing & Expansion

Suspension Culture

| | |
|---|------------------------------|
| HEKin1™ For suspension culture w/o Sodium bicarbonate, Phenol red and L-Glutamine Contains basal medium in powder form | SFM039AP-10L SFM039AP-50L |
|---|------------------------------|

Adherent Culture

| | |
|---|------------------------------|
| HEKin1™ For adherent culture w/o Sodium bicarbonate, Phenol red and L-Glutamine Contains basal medium in powder form | SFM038AP-10L SFM038AP-50L |
|---|------------------------------|

| | |
|---|------------------------|
| L-Glutamine (From non-animal source) Cell Culture Tested | TC243-10G TC243-25G |
|---|------------------------|

| | |
|--|-----------------------------|
| L-Glutamine 200mM Solution L-Glutamine in 0.85% normal saline Cell culture tested | TCL012-20ML TCL012-100ML |
|--|-----------------------------|

| | |
|--|--------------------------|
| Sodium bicarbonate powder (Sodium hydrogen carbonate) Cell culture tested | TC230-100G TC230-500G |
|--|--------------------------|

| | |
|----------------------------------|------------------------------|
| 7.5% Sodium bicarbonate solution | TCL013-100ML TCL013-500ML |
|----------------------------------|------------------------------|

Dissociation

| | |
|---|------------------------------|
| Trypsin – EDTA solution 1X 0.25% Trypsin and 0.02% EDTA in Dulbecco's Phosphate Buffered Saline w/o Phenol red | TCL007-100ML TCL007-500ML |
|---|------------------------------|

| | |
|---|--------------|
| Trypsin Inhibitor from Soyabean 1X w/ 1mg/ml of Trypsin inhibitor in Dulbecco's Phosphate Buffered Saline | TCL068-100ML |
|---|--------------|

Culture Vessels

| | Total culture area (cm ²) | Recommended volume (ml) | |
|------------------------|---------------------------------------|-------------------------|------------------------------|
| HiFactory™, 1 chamber | 647 cm ² | 200 ml | TCP204-4x1NO TCP204-8x1NO |
| HiFactory™, 2 chamber | 1279 cm ² | 400 ml | TCP205-4x1NO TCP205-8x1NO |
| HiFactory™, 5 chamber | 3175 cm ² | 1000 ml | TCP206-2x1NO TCP206-4x1NO |
| HiFactory™, 10 chamber | 6335 cm ² | 2000 ml | TCP207-2x1NO TCP207-4x1NO |

| | Surface area (cm ²) | Total Volume (ml) | |
|--|---------------------------------|-------------------|-----------------------------|
| Tissue Culture Flask Vented cap | 182 cm ² | 600 ml | TCG8-4x5NO TCG8-8x5N |
| Tissue Culture Roller Bottle Close cap | 750cm ² | 2000ml | TCG9-4x1NO TCG9-12x1NO |
| Tissue Culture Roller Bottle Vented cap | 750cm ² | 2000ml | TCG10-4x1NO TCG10-12x1NO |
| Tissue Culture Roller Bottle Close cap | 850cm ² | 2000ml | TCG15-4x1NO TCG15-12x1NO |
| Tissue Culture Roller Bottle, Expanded Surface Close cap | 1900cm ² | 2000ml | TCG16-4x1NO TCG16-12x1NO |
| Tissue Culture Roller Bottle Vented cap | 850cm ² | 2000ml | TCG17-4x1NO TCG17-12x1NO |
| Tissue Culture Roller Bottle, Expanded Surface Close cap | 4350cm ² | 5000ml | TCG18-12x1NO |

HEKin1™

HEK293 is one of the most commonly used human cell lines for the production of biotherapeutics proteins and a cell line of choice for vector-based viral vaccine production.

There are many advantages to using HEK293 cells, the major ones being –

- Easy to transfect
- Divide rapidly
- Can be utilized for both transient and stable expression
- Can be cultured in suspension or as a monolayer
- Scalable to higher volumes

And most significantly, their ability to produce fully human post-translational modifications.

The most recent example of usage of HEK293 platform for vaccine production is, COVISHIELD™, a COVID-19 vaccine containing recombinant SARS-CoV-2 spike (S) glycoprotein. This vaccine has been manufactured using genetically modified HEK293 cells.

All the HEK-based vaccine manufacturing processes require a highly efficient and productive serum free and animal component free medium at each step starting from cell banking to cell expansion and virus production in bioreactor.

HEKin1™ is a serum free and animal component free medium optimized for the growth and expansion of HEK293 cells under serum-free conditions. This is a complete media that will support growth of HEK293 cells without further supplementation. This media has been tested for its ability to support high-density cultures of HEK293 cells.



Scalability for use in shake flasks & bioreactors



Completely defined system eliminates variability



Consistent performance improves reproducibility



Decreased possibility of contamination by adventitious agents



Saves time with simplified purification and downstream processing



Supports high cell density, culture longevity and increased yield



Manufactured in GMP & ISO 9001 certified facility



HEKin1™

For Suspension Culture

Product Code : SFM039AP


Product Specifications

| Part | Name | Storage & Shipping temperature | Available Pack Sizes | |
|------|-------------------------------|--------------------------------|----------------------|-------|
| A | Basal Medium | 2-8°C | 10L | 50L |
| B | Lyophilized Growth supplement | | 10L | 5x10L |


Intended use

Intended for *in vitro* research and manufacturing processes only. Do not use for injection or infusion.

Product Performance



Scalability



Directions for Use

Preparation of complete medium

This medium does not require supplementation with Pluronic F-68®

| Ingredient | SFM039AP |
|---|--|
| Part A | 16.8g in 900ml cell culture grade water |
| Sodium bicarbonate | 3.1g sodium bicarbonate powder |
| L-Glutamine | 876 mg L-Glutamine powder |
| Add these ingredients one after the other in 900ml of cell culture grade water with constant stirring | |
| Part B preparation | <ul style="list-style-type: none">Reconstitute lyophilized growth supplement with 10ml cell culture grade water.Mix gently by pipetting up and down.Make sure that the powder is completely dissolved.Add entire content of Part B into 1L of Part A. |
| pH adjustment | Adjust the pH to 7.1 |
| Volume make up | Make up the volume to 1000ml with water. |
| Filter sterilization | Filter sterilize the complete medium immediately by filtering through a sterile membrane filter with porosity of 0.22 micron or less, using a positive pressure. |
| Antibiotic - Antimycotic solution | If required, 10ml of sterile antibiotic-antimycotic solution (A002) can be aseptically added to 1L of filter sterilized medium. |
| Storage | Store complete medium at 2 – 8°C until use |

Quality control

| Test | SFM039AP |
|---------------------------------------|--|
| Appearance | Part A: White to off-white homogenous powder Part B: White coloured lyophilized powder |
| Solubility | Clear solution at 16.8 g/L |
| pH without sodium bicarbonate | 5.90 – 6.50 |
| pH with sodium bicarbonate | 7.50 – 8.10 |
| Osmolality without sodium bicarbonate | 290 – 330mOsm/KgH ₂ O |
| Osmolality with sodium bicarbonate | 320 – 380mOsm/KgH ₂ O |
| Cultural response | The growth promotion capacity of the medium is assessed qualitatively by analyzing the cells from morphology and quantitatively by estimating cell counts. |
| Endotoxin content | NMT 1EU/ml |

HEKin1™

For Adherent Culture

Product Code : SFM038AP


Product Specifications

| Part | Name | Storage & Shipping temperature | Available Pack Sizes | |
|------|-------------------------------|--------------------------------|----------------------|-------|
| A | Basal Medium | 2-8°C | 10L | 50L |
| B | Lyophilized Growth supplement | | 10L | 5x10L |


Intended use

Intended for in vitro research and manufacturing processes only. Do not use for injection or infusion.

Product Performance



Scalability



Directions for Use

Preparation of complete medium

This medium does not require supplementation with Pluronic F-68®

| Ingredient | SFM038AP |
|--------------------|---|
| Part A | 13.3g in 900ml cell culture grade water |
| Sodium bicarbonate | 3.1g sodium bicarbonate powder |
| L-Glutamine | 876mg L-Glutamine powder |

Add these ingredients one after the other in 900ml of cell culture grade water with constant stirring

| | |
|----------------------|--|
| Part B preparation | <ul style="list-style-type: none">Reconstitute lyophilized growth supplement with 10ml cell culture grade water.Mix gently by pipetting up and down.Make sure that the powder is completely dissolved.Add entire content of Part B into 1L of Part A. |
| pH adjustment | Adjust the pH to 7.1 |
| Volume make up | Make up the volume to 1000ml with water. |
| Filter sterilization | Filter sterilize the complete medium immediately by filtering through a sterile membrane filter with porosity of 0.22 micron or less, using a positive pressure. |
| Antibiotic addition | If required, 10ml of sterile antibiotic-antimycotic solution (A002) can be aseptically added to 1L of filter sterilized medium. |
| Storage | Store complete medium at 2 – 8°C until use |

Quality control

| Test | SFM038AP |
|---------------------------------------|--|
| Appearance | Part A: White to off-white homogenous powder Part B: White coloured lyophilized powder |
| Solubility | Clear solution at 13.3 g/L |
| pH without sodium bicarbonate | 4.40 – 4.80 |
| pH with sodium bicarbonate | 7.00 – 7.60 |
| Osmolality without sodium bicarbonate | 250 – 290mOsm/KgH ₂ O |
| Osmolality with sodium bicarbonate | 340 – 380mOsm/KgH ₂ O |
| Cultural response | The growth promotion capacity of the medium is assessed qualitatively by analyzing the cells from morphology and quantitatively by estimating cell counts. |
| Endotoxin content | NMT 1EU/ml |

Procedure for Adaptation

Critical points

Cells used for adaptation should exhibit healthy morphology and have more than 90% viability. Cells should be in the mid-logarithmic phase of growth. It is necessary to subculture the cells at least thrice at each step, before going to the next step of adaptation. Subculturing should be performed when the cells are 70- 80% confluent.


In our experience, HEK293 adherent cells require gradual weaning for adaptation from serum- containing to serum-free conditions. Once adapted, the serum-free adherent culture can be directly adapted to serum-free suspension culture.

Gradual weaning

Gradual weaning is a slow procedure that involves decreasing the percentage of serum in the medium, thereby gradually adapting the cells to serum free conditions.

This procedure is applicable for adaptation of HEK293 adherent cells from serum containing conditions to serum –free conditions.

1. Subculture the cells from serum containing medium (**STEP 1**) and seed them in a 50:50 ratio of serum containing medium and SFM with a seeding density of $0.1-0.2 \times 10^6$ cells/ml.
2. Incubate at 37°C in a humidified atmosphere with 5% CO₂. Make provision for gas exchange by loosening the caps of flasks in case of closed caps or use vented caps.
3. Subculture once the cells are 70-80% confluent using Trypsin-EDTA solution.
4. Repeat steps 1 to 3 for three subcultures of each step of gradual adaptation.
5. Determine cell density and reseed the cells in 25:75 ratio of serum containing medium and SFM038AP (**STEP 2**).
6. After step 2 (25:75 serum containing medium: SFM) of adaptation, the cells cannot be directly sub-cultured in 100% serum free conditions. Complete withdrawal of serum may alter cell morphology and decrease the cell viability. Hence, it is very critical to maintain them in SFM supplemented with 1% serum (**STEP 3**), before serum free conditions.
7. When the cells reach 100% serum free step (**STEP 4**) of adaptation, subculture them till the cell density of $1.5-2.0 \times 10^6$ cells/ml is obtained within 8-10 days of culture. At this point, the cells are considered to be adapted to serum free conditions.



Note: The timelines mentioned above are based on in-house experiments for gradual adaptation. They may slightly vary depending on experimental conditions.

Direct Adaptation

This procedure is applicable for adaptation of HEK293 cells from existing serum free medium (Adherent and Suspension) to HEK293 serum-free HEKin1™ medium (Adherent and Suspension).

| | | |
|---|---|---|
| Cell harvesting from existing serum-free medium | Harvest the cells from existing serum-free medium in log phase of growth | |
| Centrifugation | Transfer the entire content to a centrifuge tube and centrifuge for 3 minutes at 1000 rpm. | |
| Cell density and Viability | Count the cells using hemocytometer, and make a note of cell count and viability. | |
| Resuspension | Carefully discard the supernatant by gentle aspiration without disturbing the pellet. Resuspend the pellet by pipetting gently with serological pipette, to get a homogenous mixture. | |
| | SFM038AP (Adherent) | SFM039AP (Suspension) |
| Reseeding | Seed with 0.2×10^6 /ml density in a new cell culture flask containing fresh HEKin1™ (SFM038AP) complete medium | Seed with 0.5×10^6 /ml density in a new cell culture flask containing fresh HEKin1™ (SFM039AP) complete medium. |
| Incubation | Incubate the cells at 37°C and 5% CO ₂ in static mode | These cells cannot be directly shifted to shaker condition. Incubate the cells at 37°C and 5% CO ₂ in static mode for 3 sub-cultures. |
| Sub-culturing | Monitor the cell health every day. Sub-culture once the cells reach 70-80% confluence by trypsinization. Neutralize the trypsin by adding equal volume of trypsin inhibitor. <i>Note: Change the medium alternate day until the cells reach 70-80% confluence. When seeded with 0.2×10^6 cells/ml, $1.5 - 2 \times 10^6$ cells/ml are obtained after subculturing</i> | Monitor the cells health under the microscope. Determine cell density and viability every day. Sub-culture when the density is double the seeding density. <i>Note: When seeded with 0.5×10^6 cells/ml, $1.5 - 2 \times 10^6$ cells/ml are obtained after subculturing</i> |
| Maintenance | Continue maintenance of these cells in static conditions for 3 passages | Continue maintenance of these cells in static conditions for total 3 passages. After 3 passages, these cells can be seeded with 0.5 million cells/ml density in a new shake flask containing fresh HEKin1™ (SFM039AP) complete medium. <i>Note: Cells should have more than 70% viability while being shifted to shaker condition.</i> |
| Incubation | — | For shaker condition, maintain the flasks at 37°C, 5% CO ₂ at 120 rpm. |
| Completion of Adaptation | Cells are considered to be fully adapted to HEKin1™ (SFM038AP) on completion of 3 passages and they can be used for further applications. | Cells are considered to be fully adapted to HEKin1™ (SFM039AP) on completion of 3 passages in shaker condition and they can be used for further applications. |

DOs and DON'Ts

- Do not use trypsin and inhibitor for SFM039AP
- Do not refrigerate cells after splitting, seed immediately.
- Vigorous pipetting will stress the cells and loss of viability.
- Do not allow cells to reach 100% confluency before sub-culture.
- Lower seeding densities may cause loss of cell viability.

Routine maintenance of adapted cells

HEK293 cells adapted to serum-free medium can be maintained in both the static culture & suspension culture.

Static Culture

For adherent culture

1. Seed adapted HEK293 cells in T-Flask at 0.1×10^6 /ml density.
2. Incubate at 37°C, 5% CO₂.
3. Monitor the cell health every day. Subculture once the cells reach 70-80% confluence by trypsinization. Neutralize the trypsin by adding equal volume of trypsin inhibitor.
4. Centrifuge the cells at 1000rpm. Discard the spent medium.
5. Re-suspend the pellet in fresh medium & re-seed in new vessel at 0.1×10^6 /ml density.

Shaker Culture

For suspension culture

1. Seed the adapted cells into shaker flask of required volume with seeding density of 0.5×10^6 /ml.
2. Incubate at 37°C, 5% CO₂ & 120 rpm.
3. Determine the cell density every day & subculture them when the density reaches $1.5-2.0 \times 10^6$ cells/ml.
4. Centrifuge the cells at 1000rpm. Discard the spent medium.
5. Re-suspend the pellet in fresh medium & re-seed in new vessel at 0.5×10^6 /ml density.

Bioreactor Cultivation

Users are recommended to optimize incubation density, DO, agitation & other parameters empirically, depending on bioreactor scale. Supplementation with Pluronic F68 is not required in bioreactor.

Storage and shelf life

Store basal medium (Part A) at 2-8°C away from bright light.

Store serum free lyophilized growth supplement (Part B) at 2-8°C.

Use before expiry date given on the product label.

Shelf life of the complete medium is 4-6 weeks at 2-8°C.

Note: Freezing of the basal medium and complete medium is not recommended.



For Life is Precious

HiMediaLaboratories™
www.himedialabs.com

CORPORATE OFFICE

HiMedia Laboratories Pvt. Ltd., Plot No. C40, Road No. 21Y, MIDC, Wagle Industrial Area, Thane (West) 400604, Maharashtra, India. Tel No.: 022-69034800.
Email: atc@himedialabs.com. Website: www.himedialabs.com.

