
HiViral™ Transport Medium

With Viral Transport Medium

Recommended for collection and transport of viruses

Product Code: AL167

Intended Use: For collection and transportation of samples containing Viruses Chlamydiae, Mycoplasma and Ureaplasma from the collection site to the laboratory

Product Description :

HiViral™ Transport medium is a specially formulated medium for collection, transport and long term freeze storage of viruses. It is designed to maintain the optimum viability and virulence of the viral sample.

HiViral™ Transport Medium is made of Hanks Balanced Salt Solution and contains a protective protein antibiotics to control microbial contamination and buffers to control the pH. Phenol red is used as a pH indicator. The medium also contains a cryoprotectant which helps in preserving the viruses, if specimens are frozen for prolonged storage.

Composition :

Proprietary

Quality Control:

Appearance

Orange colored, clear solution.

pH at 25°C

7.3 ± 0.3

Osmolality in mOsm/Kg H₂O

500mOsm* - 600mOsm*

Sterility

No bacterial or fungal growth was observed after 14 days of incubation as per USP specification.

Precautions:

1. Isolation of viruses will largely depend on proper specimen collection, timing of sample collection and processing of samples.
2. Specimen collection should be done in the acute phase of illness.
3. Do not use the product if,
 - (i) there is change in the color of the medium,
 - (ii) there is evidence of leakage and
 - (iii) there are other signs of deterioration.
4. To maintain infectivity of viruses, it is important that temperature be properly maintained for sample collection to processing.
5. Avoid repeated freeze-thaw of collected samples.
6. It is recommended to refer to standard procedures and published protocols for sample collection and processing.

Storage and Shelf Life:

HiViral™ Transport medium should be stored at 15-30°C.

Use before expiry date given on the product label.

IVD

In vitro diagnostic medical device



CE Marking



Consult instructions for use



Do not use if package is damaged



Reg. Off : 23, Vadhani Ind Est.,
LBS Marg, Mumbai-400086, India.
Works : B-4-5-6 / MIDC, Palkhed,
Dindori, Nashik- 422202
Maharashtra, India
www.himedialabs.com

EC REP

CE Partner 4U, Esdoornlaan 13, 3951
DB Maarn The Netherlands,
www.cepartner4u.eu



Single use.
Not intended to be
reprocessed and/or used on
another patient

Revision : 5 / 2020

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.