



HiViral™ Nasopharyngeal Nylon Flocked Swab

With breakpoint

Product Code: PW1172

Specifications

Tip

Material of construction	: Nylon fiber
Flocked	: Yes
Diameter	: 3.4mm
Length	: 18.4mm
Length of fibre	: 1.04mm

Shaft

Length	: 153mm
Breakpoint	: Yes
Breakpoint at	: 89mm from top

Sterilization	: Gamma irradiation
Packing	: Individually packed in Glassene poly pouch

Features and Applications

The ultra-flexible plastic shaft is ergonomically designed to optimize the collection efficiency with improved patient comfort. The capillary action between the nylon fiber strands facilitates strong hydraulic uptake of liquid samples. The sample stays close to the surface allowing easier and rapid elution, which results in minimum hands-on specimen processing.

- Suitable for nasal and nasopharyngeal sampling
- Available in unique molded breakpoint design