

TC019 Amphotericin B, Cell Culture Tested**Product Number****Packing**

TC019	: 1G
TC019	: 5G

Product Information

Product Code	: TC019
Product Name	: Amphotericin B, Cell Culture Tested
Synonym	: Fungizone
Molecular Formula	: C ₄₇ H ₇₃ NO ₁₇
Molecular Weight	: 924.08
CAS No.	: 1397-89-3
EC No.	: 215-742-2
HS Code	: 2941 90 00
Shelf Life	: 3 years


Technical Specification

Appearance	: Yellow to orange crystals or powder
Solubility	: Clear yellow solution at 3mg in 1ml of dimethyl sulphoxide
Cell Culture Test	: Passes
Antibiotic sensitivity	: Passes
Potency	: min. 750 µg/mg
Assay (HPLC)	: min. 90.00%

GHS Safety Information

Hazard Statement(s)	: H315-H319-H335
Precautionary Statement(s)	: P261-P305+P351+P338
Signal Word	: Warning
Hazard Pictogram(s)	

Risk and Safety Information

R-Phrase(s)	: 36/37/38
S-Phrase(s)	: 26-36
WGK	: 3
RTECS	: BU2625000
Storage Temperature(°C)	: Store at 2 - 8°C
Hazard Symbol(s)	

Transport Information

Marine Pollutant	: No
------------------	------



Product Information

ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such Information or with respect to any product, method or apparatus referred to herein