

**TC202 Vancomycin hydrochloride, Cell Culture Tested****Product Number Packing**

TC202	: 1G
TC202	: 500MG
TC202	: 5G

**Product Information**

Product Code	: TC202
Product Name	: Vancomycin hydrochloride, Cell Culture Tested
Molecular Formula	: $C_{66}H_{75}Cl_2N_9O_{24}.HCl$
Molecular Weight	: 1485.71
CAS No.	: 1404-93-9
HS Code	: 2941 90 00
Shelf Life	: 3 years

**Technical Specification**

Appearance	: White to brown powder
Solubility	: Clear faint yellow to faint orange solution at 50mg in 1ml of water
pH (5% in water at 25°C)	: 2.50 - 4.50
Antibiotic sensitivity	: Complies
Cell Culture Test	: Passes
Potency	: min. 900 µg/mg
Assay (HPLC)	: min. 80.00%

**GHS Safety Information**

Hazard Statement(s)	: H317
Precautionary Statement(s)	: P280
Signal Word	: Warning
Hazard Pictogram(s)	

**Risk and Safety Information**

R-Phrase(s)	: 43
S-Phrase(s)	: 36/37
WGK	: 2
RTECS	: YW4380000
Storage Temperature(°C)	: Store at 2 - 8°C
Hazard Symbol(s)	

**Transport Information**

Marine Pollutant	: No
------------------	------



## Product Information

ADR/RID : Not Dangerous Goods  
IMDG : Not Dangerous Goods  
IATA : Not Dangerous Goods

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919  
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

---

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such Information or with respect to any product, method or apparatus referred to herein