

---

## Tylosin tartrate, Cell Culture Tested

**TC205**

---

### Product Identifier

EC No.	:	215-781-5
HS Code	:	2941 90 19
Storage	:	2 - 8°C
Shelf life	:	3 years


---

### Technical Specification

Appearance	:	Off-white to yellow hygroscopic powder
Solubility	:	33.3 mg soluble in 1 mL of water
pH (2.5 % in water at 25°C)	:	5.00 - 7.20
Potency	:	min. 900 µg/mg
Cell Culture Test	:	Passes

---

### Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H317- H334
Precautionary Statement(s)	:	P261- P280- P285- P342 + P311- P363- P501
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	3