

Krebs-Ringer Bicarbonate Buffer 1X

**With 1.8gms Glucose per liter and Sodium bicarbonate
Without Calcium Chloride**

Product Code: TL1097

Product Description :

All media used in tissue culture have a basis of a synthetic mixture of inorganic salts known as a physiological or balanced salt solution (BSS). All the physiological salt solutions have been derived from the salt solution originally described by Sydney Ringer (1885). The first balanced salt solution to be developed specifically for supporting the metabolism of mammalian cells was Tyrode's solution. Since then many modifications have been done to obtain better buffering salt solutions and to prevent calcium precipitation.

The function of salt solution is:

- To maintain the medium within physiological pH range.
- To maintain intracellular and extra cellular osmotic balance.
- Modified with a carbohydrate, such as glucose serves as an energy source for cell metabolism.

TL1097 is Krebs-Ringer Bicarbonate Buffer with 1.8gms Glucose per liter and sodium bicarbonate.

Composition :

Ingredients	mg/L
INORGANIC SALTS	
Disodium hydrogen phosphate anhydrous	100.000
Magnesium chloride anhydrous	46.800
Potassium chloride	340.000
Sodium bicarbonate	1260.000
Sodium chloride	7000.000
Sodium dihydrogen phosphate anhydrous	180.000
OTHERS	
D-Glucose	1800.000

Revision: 02/2019

Quality Control:

Appearance

Clear, colourless solution

pH

7.00 -7.60

Osmolality in mOsm/Kg H₂O

240.00 -280.00

Sterility

No bacterial or fungal growth is observed after 14 days of incubation, as per USP specification.

Toxicity test

Passes

Endotoxin Content

NMT 1 EU/ml

Storage and Shelf Life:

Store at 15-30°C away from bright light.

Shelf life is 24 months.

Use before expiry date given on the product label.

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ Publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.